




HtrA4 ELISA

REF Cat. No.: 30 516 401

For quantitative determination of human HtrA4

 96 Determinations
 2°C – 8 °C



For research use only
Not for diagnostic use



Contents

Short Review	3
Assay description	4
Principle of the assay	5
Safety warning and precautions	6
Storage	6
Components of the assay system	6
Sample preparation and storage	7
Critical parameters	7
Preparation of reagents	8
Assay protocol	9
Protocol summary	10
Data processing	11
Additional information	12
Specificity	12
Sensitivity	12
Linearity	12
Recovery	12
Reproducibility	12
Troubleshooting	13
Related Products	14
Distributors	15



BioTeZ Berlin-Buch GmbH
Robert-Rössle-Str. 10, Haus 72
13125 Berlin
Tel.: ++49 (0)30-9489 3322
Fax: ++49 (0)30-949 2008
e-mail: pro@biotez.de
<http://www.biotez.de>

This manual is valid from April, 2016








Short Review

Human HtrA4, the newest member of the HtrA family, an unique family of multidomain serine proteases, was initially reported as a serine protease associated with pregnancy [1]. Polypeptide chains of human HtrA4 consist of several domains: An N-terminal insulin-like growth factor domain, followed by a Kazal-type serine protease inhibitor domain, a linker region, a trypsin-like protease domain and a PDZ domain [2].

Two studies reported in 2012 an association between abnormal levels of HtrA4 and preeclampsia [1, 3]. However, the findings are contradictory, while the study by Inagaki *et al.* [3] reported an elevation of HtrA4 in preeclampsia, the study by Wang *et al.* [1] published that HtrA4 expression was decreased in preeclampsia. Therefore, the precise association between HtrA4 and preeclampsia was for some years unclear, till recently, when Sing *et al.* [4] published new data that showed that human HtrA4 expression was restricted to the placenta and significantly up-regulated in early-onset but not late-onset preeclampsia. The serum HtrA4 levels in normal pregnancy increased significantly between the first and second trimesters, then remained constant. Women with early-onset but not late-onset preeclampsia showed significantly higher HtrA4 levels in serum compared with gestational age-matched controls [4].

References

1. Wang LJ, Cheong ML, Lee YS, Lee MT, Chen H. High-temperature requirement protein A4 (HtrA4) suppresses the fusogenic activity of syncytin-1 and promotes trophoblast invasion. *Mol Cell Biol.* 2012; 32(18):3707–3717
2. Zumbunn, J. and Trueb, B. (1996) FEBS Lett. 398, 187-192
3. Inagaki A, Nishizawa H, Ota S, Suzuki M, Inuzuka H, Miyamura H, Sekiya T, Kurahashi H, Udagawa Y. Upregulation of HtrA4 in the placentas of patients with severe pre-eclampsia. *Placenta.* 2012; 33(11):919–926.
4. Singh, H., Zhao, M., Chen, Q., Wang, Y., Li, Y., Kaitu'u-Lino, T., Tiong, S., Nie, G., Human HtrA4 Expression Is Restricted to the Placenta, Is Significantly Up-Regulated in Early-Onset Preeclampsia, and High Levels of HtrA4 Cause Endothelial Dysfunction, *JCEM* 100, 7, (2015), doi: 10.1210/jc.2014-3969

REF	Catalogue-No.	LOT	Batch-No.
	Storage temperature		Producer
	Notice advices		Use by
	Consult Instructions for use		Number of determination
			Biohazard

Assay Description


The HtrA4 ELISA from BioTeZ Berlin-Buch GmbH provides a highly sensitive and specific quantitative determination of HtrA4 in serum, tissue (Placenta) and cell culture supernatants.

The calibration curve covers the range from 80 ng/ml to 1.25 ng/ml.

This quantitative assay is based on a two site sandwich format. A highly specific monoclonal antibody against HtrA4 is immobilised on the plate. HtrA4 will be bound to the wells, other components of the sample are removed by discarding/drying by tapping and washing of the plate. The analyte is detected in two steps using a secondary biotin-labeled monoclonal antibody and a highly polymerised streptavidin-peroxidase conjugate. Any excess is removed by discarding/drying by tapping and washing after each detection step. The amount of peroxidase bound to each well is determined by the addition of TMB Substrate. The reaction is stopped by adding the Stop Solution and the yellow colour is read in a microtiter plate reader at 450 nm. The concentration of HtrA4 in a sample is determined by interpolation from the standard curve.

One kit contains reagents for 96 determinations, thus allowing the measurement of one standard curve and 40 samples in duplicate.

Principle of the assay



Step 1: Incubation of standard or sample on the microtiter plate. Specific binding of HtrA4
(duration: 120 minutes)

Step 2: Detection of bound HtrA4 with biotinylated antibody
(duration: 90 minutes)

Step 3: Addition of streptavidin-peroxidase conjugate
(duration: 30 minutes)

Step 4: Color development after addition of TMB Substrate
(duration: 5-15 minutes)

Fig. 1: Scheme of the assay procedure

Safety warning and precautions

Warning: Stop Solution contains 0.25 M sulfuric acid. Wear eye, hand, face and clothing protection when using this material.

The Serum Standard Diluent contains human serum that has been negatively tested for contaminations of HAV, HBV, HCV and HIV. Wear suitable protective clothing, such as laboratory overalls and gloves and observe caution when working with this material.

All chemicals should be considered as being potentially hazardous. We therefore recommend that this product is handled only by those persons who have been trained in laboratory techniques and that it is used in accordance with the principles of good laboratory practice. Wear suitable protective clothing, such as laboratory overalls, safety glasses and gloves. Avoid contact with skin and eyes. In case of skin and eyes contact, wash immediately with water.

Storage

All components of the kit can be stored in the refrigerator (2-8 °C).

Note! Wash buffer concentrate can precipitate in refrigerator. We recommend stored at 18-35°C

HtrA4 Standard Puffer should be used immediately or stored at -20 °C.

Once reconstituted, the **HtrA4 Standard** should be used immediately or stored at -80 °C.

Diluted **Biotinylated antibody** and diluted **Conjugate solution** should be prepared freshly directly before use. When running a partial plate, only suitable aliquots of these solutions should be made.

Components of the assay system

Microtiter plate: The plate contains 6 x 16 strips coated with monoclonal anti-HtrA4.

Assay Buffer 2.5: Bottle contains 25 ml buffer with additives.

Wash Buffer Concentrate: Two bottles containing each 25 ml of a buffer concentrate. It has to be diluted 20-fold with distilled or deionised water before use.

HtrA4 Standard (lyophilised): The vial contains 400 ng of lyophilized HtrA4. It has to be reconstituted with 1 ml of HtrA4 Standard Buffer prior to use to get a 400 ng/ml standard solution.

HtrA4 Standard Buffer: Solution to reconstitute HtrA4 lyophilized.

Serum Standard Diluent: This solution (20 ml) contains certified human serum in a buffer with additives. Ready to use.

Warning: Component contains human serum! Wear eye, hand, face and clothing protection when using this material!

Detection Buffer: Solution (25 ml) consists of phosphate buffer containing additives.

Biotinylated Antibody: Solution (150 µl) contains biotinylated monoclonal anti-HtrA4 antibody in buffer with additives. The antibody has to be diluted 100-fold with Detection Buffer before use.

Conjugate Solution: This solution (300 µl) contains a highly polymeric streptavidin-peroxidase conjugate with preservatives. It has to be diluted 100-fold with Detection Buffer prior to use.

TMB Substrate: Bottle contains 12 ml of a TMB solution ready to use.

Stop Solution: Bottle contains 12 ml 0.25 M H₂SO₄. Ready to use. **Warning:** Stop Solution contains 0.25 M sulfuric acid. Wear eye, hand, face and clothing protection when using this material!

EQUIPMENTS REQUIRED BUT NOT PROVIDED

- Pipettes with disposable tips (100 µl, 500 µl and 1 ml), a multi-channel pipette (100 - 200 µl) would be appropriate
- Distilled or deionised water
- Horizontal orbital microplate shaker
- Microplate reader capable of measuring at 450 nm

Sample preparation and storage

Serum:

- Serum samples may be stored at -80°C. When stored at -80°C, it is absolutely necessary to mix the samples thoroughly prior to measuring. Avoid freeze-thaw cycles.
- Dilute the serum samples minimum **1:10** with Serum Standard Diluent, depending on the possible concentration of the analyte.

Cell culture supernatant

- Centrifuge the samples to remove any particles. The supernatants can be stored at -80°C. Avoid freeze-thaw cycles.
- Dilute the samples minimum **1:10 or more with Assay Buffer 2.5**, depending on the possible concentration of the analyte. For measuring these samples, use the **standards prepared with Assay Buffer 2.5**.

Critical Parameters

- Allow samples and all reagents to equilibrate to room temperature (18-35°C) prior to performing the assay. This is especially a prerequisite for Wash Buffer and TMB Substrate!
- For the highly sensitive determination of HtrA4 in serum, use the Serum Standard Diluent to prepare your standard curve. For all other samples, use Assay Buffer 2.5.
- It is absolutely important that all wells are washed thoroughly and uniformly. When washing is done by hand, ensure that all wells are completely filled and emptied at each step by discarding the content of the plate with forceful motion and drying by tapping the inverted plate on dry absorbent surface (see technical advice on the following link: <http://www.youtube.com/watch?v=FZirnCas17Y>)
- Use only reagents from the same lot for each assay. This is especially important when running more than one plate per sample group.
- A separate standard curve must be run on each plate.
- Mix all reagents thoroughly prior to use, but **avoid foaming!**
- Keep the wells sealed with the foil except when adding reagents and during reading.
- Any variation in the protocol can cause variation in binding!
- The kit should not be used beyond the expiration date on the kit label.
- The values obtained by the samples should be within the standard range. If this is not the case, dilute the sample and repeat the assay.
- We take great care to ensure that this product is suitable for all validated sample types, as designated in this manual. However, it is possible that in some cases, high levels of interfering substances may cause unusual results.
- EDTA (0.5 mM - 4 mM) **does not** inhibit the assay!

Preparation of reagents

Please note: To prevent margin effects it is absolutely necessary to equilibrate all reagents **to room temperature (1 h incubation on the bench)** prior to use. For the dilution of the wash buffer concentrate use either distilled or deionised water. Always seal the plates with the provided foil during incubation!

Assay Buffer 2.5: After equilibration to room temperature, the buffer is ready to use.

Wash Buffer: Dilute the contents of a 25 ml bottle to 500 ml with distilled or deionised water. Make sure that the whole content of the bottle is used by repeated washing.

HtrA4 Standard (lyophilised): Add 1 ml HtrA4 Standard Buffer to the standard tube (yellow lid) and allow the contents to dissolve for 5-10 minutes. Gently mix, but avoid foaming of the reagent!

HtrA4 Standard Buffer: After equilibration to room temperature, the buffer is ready to use.

Serum Standard Diluent: After equilibration to room temperature, the reagent is ready to use.

Detection Buffer: After equilibration to room temperature, the reagent is ready to use.

Biotinylated Antibody: Dilute Biotinylated Antibody 100-fold with Detection Buffer. For a whole plate add 120 μ l from the antibody solution tube (red lid) to 12 ml Detection Buffer. When running half a plate, add 60 μ l antibody solutions to 6 ml Detection Buffer.

Conjugate Solution: Dilute the provided Conjugate Solution 100-fold with Detection Buffer.

For the whole plate, add 120 μ l of the conjugate tube (blue lid) to 12 ml Detection Buffer. When running half a plate, add 60 μ l of the conjugate to 6 ml of Detection Buffer.

TMB Substrate and Stop Solution: After equilibration to room temperature, the reagents are ready to use.

Preparation of standards with Assay Buffer 2.5

(For cell culture supernatants)

1. Label 7 tubes with 1.25, 2.5, 5, 10, 20, 40 and 80 ng/ml.
2. Pipette 400 μ l of Assay Buffer 2.5 into the 80 ng/ml tube, in the remaining tubes pipette 250 μ l of Assay Buffer 2.5.
3. Pipette 100 μ l of the stock standard (400 ng/ml) into the 80 ng/ml tube and mix thoroughly.
4. Pipette 250 μ l of the 80 ng/ml standard into the tube labelled with 40 ng/ml and mix thoroughly.
5. Repeat this dilution procedure with the other standard tubes.
6. The blank value (0 ng/ml) is obtained by using only Assay Buffer 2.5.
7. The stock solution is not part of the standard curve and can be stored at -20°C

Preparation of standards with Serum Standard Diluent

(For the highly sensitive measurement of serum samples)

1. Label 7 tubes with 1.25, 2.5, 5, 10, 20, 40 and 80 ng/ml.
2. Pipette 400 μ l of Serum Standard Diluent into the 80 ng/ml tube, in the remaining tubes pipette 250 μ l of Serum Standard Diluent.
3. Pipette 100 μ l of the HtrA4 standard (400 ng/ml) into the 80 ng/ml tube and mix thoroughly.
4. Pipette 250 μ l of the 80 ng/ml standard into the tube labelled with 40 ng/ml and mix thoroughly.
5. Repeat this dilution procedure with the other standard tubes.
6. The blank value (0 ng/ml) is obtained by using only Serum Standard Diluent.
7. The stock solution is not part of the standard curve and can be stored at -20°C.


Assay protocol

1. Prepare reagents and standards as described in the sections above. For assessing serum samples, prepare the standard using the Serum Standard Diluent. **Remind that it is necessary to equilibrate the reagents to room temperature before use.**
2. Prepare the unknown samples as described above by appropriate dilution with Assay Buffer 2.5.
3. Prepare the Microtiter plate by inserting the required amount of wells into the frame. **Note that you need 16 wells for the standard curve.**
4. Pipette 100 μ l of the reconstituted standards (80 ng/ml) in duplicate in the wells using a clean pipette tip for each standard. Assay Buffer 2.5 serves as zero blank.
5. Pipette 100 μ l of the prepared unknown samples in duplicate into the wells.
6. Seal the plate with the provided foil and incubate on a shaker at room temperature for exactly 120 minutes.
7. Wash by filling each well with Wash Buffer (200 μ l), then remove by discarding/drying by tapping inverted plate against clean paper towels. Take care that all wells are completely filled and emptied at each wash. Wash the wells 4 times with Wash Buffer.
8. Add 100 μ l of diluted Biotinylated Antibody into each well.
9. Seal the plate and incubate on a shaker at room temperature for exactly 90 minutes.
10. Wash by filling each well with Wash Buffer (200 μ l), then remove by discarding/drying by tapping inverted plate against clean paper towels. Take care that all wells are completely filled and emptied at each wash. Wash the wells 4 times with Wash Buffer.
11. Add 100 μ l of diluted Conjugate Solution into each well.
12. Seal the plate with the provided foil and incubate on a shaker at room temperature for exactly 30 minutes.
13. Wash by filling each well with Wash Buffer (200 μ l), then remove by discarding/drying by tapping inverted plate against clean paper towels. Take care that all wells are completely filled and emptied at each wash. Wash the wells 6 times with Wash Buffer.
14. Add 100 μ l of TMB Solution to each well.
15. Seal the plate with foil provided and incubate in the dark at room temperature without shaking for 10 minutes.
16. Stop the reaction by adding 100 μ l of Stop Solution to each well.
17. Read the plate at 450 nm (620 nm reference filter) within 30 minutes. Reading of the plate without reference may yield higher absorbance and thus may be less accurate.


Protocol summary

Prepare reagents, standards and samples as described, equilibrate reagents to **room temperature**.


For assessing human serum samples, prepare a Serum Standard!




Pipette 100 μ l standard or sample in duplicate into the wells. Incubate 120 minutes at room temperature on a shaker




Discard/Dry by tapping/Wash 4 times




Add 100 μ l of diluted Biotinylated Antibody to each well. Incubate 90 minutes at room temperature on a shaker.




Discard/Dry by tapping/Wash 4 times.




Add 100 μ l of diluted Conjugate Solution to each well. Incubate for 30 minutes at room temperature on a shaker.



Discard/Dry by tapping/Wash 6 times.



Add 100 μ l of TMB Substrate to each well. Incubate for 10 minutes at room temperature **in the dark** without shaking!



Add 100 μ l Stop Solution to each well. Read at 450 nm within 30 minutes (when possible with reference filter at 620 nm).

Data processing

Calculation of results

The calculation is illustrated using representative data: the assay data should be similar to that shown in table 1.

1. Calculate the average absorbance for each set of standard wells.
2. A standard curve is generated by plotting the mean absorbance (x-axis, fig. 2) against ng/ml standard (y-axis, fig.2).
3. The ng/ml values of the samples can be read directly from the graph or calculated by the regression coefficients.
4. Multiply the calculated ng/ml values by the dilution factor of the samples.

Standard (ng/ml)	Standard curve in Assay Buffer 2.5	
	Absorbance (450 nm)	Average
80.00	2.343 / 2.474	2.409
40.00	1.359 / 1.337	1.348
20.00	0.698 / 0.653	0.676
10.00	0.319 / 0.340	0.330
5.00	0.184 / 0.184	0.185
2.50	0.120 / 0.118	0.120
1.25	0.087 / 0.083	0.085
0	0.052 / 0.050	0.051

Table 1: Typical assay data




Fig. 2: Typical standard curve (prepared with Assay Buffer 2.5).

Standard (ng/ml)	Standard curve in Serum Standard Diluent	
	Absorbance (450 nm)	Average
80.00	2.095 / 2.227	2.161
40.00	1.174 / 1.124	1.150
20.00	0.572 / 0.561	0.567
10.00	0.278 / 0.280	0.280
5.00	0.154 / 0.155	0.155
2.50	0.105 / 0.105	0.104
1.25	0.073 / 0.072	0.073
0	0.046 / 0.047	0.047

Table 2: Typical assay data




Fig. 3: Typical standard curve (prepared with Serum Standard Diluent).

Additional Information

Specificity

The HtrA4 ELISA has a high sensitivity and high specificity for quantitative determination of HtrA4. The assay recognizes HtrA4 full length (50 kDa) and the truncated form of HtrA4 (138-476 aa, 37 kDa). **It does not cross-react with human HtrA1 or HtrA3.**

Sensitivity

The sensitivity was determined up to 1 ng/ml.

Linearity

The following samples below were measured after dilution with Assay to assess linearity of the assay.

Sample Type	Added concentration (ng/ml)	Measured (ng/ml)	Recovery (%)
Assay Buffer	80	79.8	100
	40	40.4	101
	20	19.7	98
	10	9.3	93

Recovery

The recovery of HtrA4 standard spiked to levels throughout the range of the assay in serum.

Matrix	Recovery range (%)	Average (%)
Serum (n=6)	84-104%	94%

Reproducibility

Within assay precision

The within assay precision was measured by assaying four control samples in 10 duplicates in Assay Buffer on one plate

Standard (ng/ml)	Standard deviation	%CV	N
80	0.07	3.36	10
40	0.07	5.22	10
20	0.02	3.44	10
10	0.01	3.84	10

Between assay precision

The between assay variation was measured by assaying

Standard (ng/ml)	Standard deviation	%CV	N
80	0.11	4.3	10
40	0.08	5.3	10
20	0.11	4.5	10
10	0.03	7.2	10

Troubleshooting

Problem	Potential cause	Recommendation
Low absorbance	<ul style="list-style-type: none">• Wrong wavelength• Enzyme conjugate out of date/reagents improperly stored• Improper incubation time and temperature• Reagents not equilibrated to RT• Reagents not correctly prepared	<ul style="list-style-type: none">• Check reader wavelength• Control the expiration date/storage conditions• Control the incubation time and temperature• Check equilibration of reagents to RT• Check preparation of reagents
High absorbance/ high zero standard value (>0.15 OD)	<ul style="list-style-type: none">• Incomplete washing• Improper removing of residual fluid• Improper incubation time and temperature• Reagents not equilibrated to RT• Reagents not correctly prepared	<ul style="list-style-type: none">• Ensure that every well is completely filled/emptied during each washing step• Check that plates are blotted on tissue paper after each washing step• Control the incubation time and temperature• Check equilibration of reagents to RT• Check preparation of reagents
Flat curve/poor reproducibility	<ul style="list-style-type: none">• Wrong wavelength• Enzyme conjugate out of date/reagents improperly stored• Improper preparation of working standards• Pipette errors• Contamination of components by use of unclean reservoirs/used pipette tips• Margin effects by using of cold substrate solution• Washing incomplete	<ul style="list-style-type: none">• Check wavelength• Control the expiration date/storage conditions• Check preparation of standards• Check pipette calibration• Use separate reservoirs and always new pipette tips• Equilibrate substrate to room temperature• Ensure sufficient washing procedure

Related Products

For enrichment of proteases you can use:

Cat.No.: 30 501 101 EXTRACT HtrA1 mini
Cat.No.: 30 501 102 EXTRACT HtrA1 maxi
Cat.No.: 30 501 501 EXTRACT HtrA3 mini
Cat.No.: 30 501 502 EXTRACT HtrA3 maxi
Cat.No.: 30 501 401 EXTRACT HtrA4 mini
Cat.No.: 30 501 402 EXTRACT HtrA4 maxi

For detection of proteases you can use our antibodies, like:

Cat.No.: 31 601 004 anti-HtrA1
Cat.No.: 31 605 004 anti-HtrA3
Cat.No.: 31 604 004 anti-HtrA4

More antibodies information you can find on our catalogue:

(http://www.biotez.de/images/biotez/downloads/Catalogue_BioTeZ_2016.pdf)

For determination of concentration, you can use our ELISA kits:

Cat.No.: 30 516 101 HtrA1 ELISA
Cat.No.: 30 516 501 HtrA3 ELISA
Cat.No.: 30 518 101 MMP13 ELISA

For determination of protease activity, you can use our activity ELISAs:

Cat.No.: 30 510 111 Activity Aggrecanase Assay
Cat.No.: 30 510 211 Activity Sensitive Aggrecanase Assay
Cat.No.: 30 516 103 Universal Protease Activity ELISA

For more products, please use our catalogue.

Distributors

United States

Eagle Biosciences, Inc
Dan Keefe
10 Columbia Drive
AMHERST, NH 03031
Tel: +617-419-2019
Fax: +617-419-1110
<http://www.eaglebio.com>

United States

B-Bridge International, Inc.
320 Logue Ave
Mountain View
CA 94943 USA
Tel: +650-969-7727
Fax: +650-969-7737

Italy

Divisione BIOSPA della SPA
Società Prodotti Antibiotici S.p.A.
Via Biella 8 - 20143 Milano – ingresso
Via Modica 6
Tel: + 39-02-891 39 545
Fax: + 39-02-891 20 996

Japan

B-Bridge Co., Ltd.
Yashima Building 5F
1-4-5 Sarugaku-cho,
Chiyoda-ku, Tokyo 101-0064 Japan
Tel: +81-3-5281-4611
Fax: +81-3-5281-4612
<http://lifescience.b-bridge.com>



BioTeZ Berlin-Buch GmbH

Robert-Rössle-Str. 10
13125 Berlin Germany
Tel.: ++49 (0)30-9489 3317
Fax: ++49 (0)30-949 4509
e-mail: pro@biotez.de
<http://www.biotez.de>